careCL TECH FIT

User Guide / Guía del usuario

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1 Product Introduction

AGELESS BEAUTY

careCL originated from the **experience with dermatology aesthetic devices** and the expertise of professionals.

To provide an alternative to the high cost of procedures and the inconvenience of frequent clinic visits, it was developed so that safe and convenient skin firmness care can be performed at home.

careCL **TECH FIT** is a home beauty device designed based on energy-based technology used in professional aesthetic care, allowing users to manage skin firmness conveniently at home.

Rather than simply delivering heat, the device is designed with a structure that delivers controlled energy to the skin safely and evenly.

2 Package Components

careCL **TECH FIT**
device

Grid Mask Pack
(3 sheets)

USB-C
Charging Cable

User Guide &
Warranty Card

careCL
Collagen Gel

Cradle

3

Product Overview

careCL **TECH FIT** consists of three care modes: SHOT / GLOW / EMS, which can be selected depending on skin condition and care purpose.

Category	SHOT	GLOW	EMS
Purpose	Uses stamping energy-based technology to target specific skin areas and helps support skin firmness and improve the appearance of wrinkles.	Helps improve skin radiance and smoothness, promoting a healthy-looking skin tone.	Uses micro-current to gently stimulate facial muscles, helping to support facial contours, skin firmness, and a refreshed appearance.
Treatment Areas	Face, eye area, neck, arms, abdomen	Entire face	Facial contour, neck
Intensity Level	Level 1-5	Level 1-5	Level 1-5
Vibration Function	Basic function	ON/OFF	ON/OFF
Recommended Frequency	2-3 times per week	Daily use possible	Daily use possible
Operating Time	7 minutes	5 minutes	5 minutes
Safety Features	Smart contact sensor (operates only when in contact with skin), auto stop, auto power OFF	Contact sensor, auto stop, auto power OFF	Contact sensor, auto stop, auto power OFF
User Guide	When using with Grid Mask Pack, starting from LEVEL 3 is recommended. When using with Collagen Gel, starting from LEVEL 1 is recommended.	For sensitive skin, starting from LEVEL 1 is recommended.	For sensitive skin, starting from LEVEL 1 is recommended.
LED	FACE SHOT RED supports skin firmness care and improvement in the appearance of wrinkles EYE SHOT GREEN supports a brighter and more even-looking skin tone BODY SHOT YELLOW supports skin comfort care and helps reduce the appearance of puffiness	PURPLE supports overall skin condition and skin vitality	RED PURPLE YELLOW GREEN Auto Cycle




4 Treatment Areas

* This product is not a medical device. Do not use for medical treatment purposes.

Mode 1 SHOT 5-Second Stamping Shot


SHOT mode is designed with a "5-second stamping shot" technique that focuses energy on the skin by maintaining contact with one area, rather than rubbing across the skin.

careCL's Patented Technology
(5-Second Stamping Shot)



* Do not use on the following areas: Eyeballs, eyelids, lips, uvula area, thyroid area

Mode 2 GLOW

Mode 3 EMS

Mode Selection Guide

Basic Care Routine

SHOT mode (Skin Firmness Care) ➔ **GLOW mode** (Skin Texture & Skincare Absorption Care)
➔ **EMS mode** (Skin Care Finishing)

Recommended order : **Skin Firmness Care** ➔ **Skin Texture Care** ➔ **Finishing Care**

Total Care Time: Approximately 10–20 minutes.

* Adjust the treatment time depending on your care needs.

MODE	ENERGY SOURCE	TIME	FREQUENCY	HOW TO USE / EXPECTED BENEFITS
SHOT	FACE Energy-based technology (3,2,1Mhz) RED LED , 3D Vibration	7 minutes	3 times per week	FACE SHOT Mode Apply skincare product to the skin, then place the device head on the target area and hold for stamping (5 seconds). Helps support skin firmness and smoother-looking skin.
	EYE Energy-based technology (3,2Mhz) GREEN LED , 3D Vibration	7 minutes	3 times per week	EYE SHOT Mode Apply skincare product to the skin, then place the device head on the eye area and hold for stamping (5 seconds). Supports smoother-looking skin around the eye area and a brighter-looking skin tone.
	BODY Energy-based technology (2,1Mhz) YELLOW LED , 3D Vibration	7 minutes	3 times per week	BODY SHOT Mode Apply skincare product to the skin, then place the device head on the target body area and hold for stamping (5 seconds). Supports skin care for firmer-looking skin.
GLOW	Galvanic Ion Electroporation VIOLET LED , 3D Vibration (Optional)	5 minutes	Daily	GLOW Mode Apply skincare product to the skin and gently glide the device across the treatment area. Helps improve skin radiance, support skincare absorption, and promote smoother-looking skin and refined-looking pores.
EMS	EMS(Low Frequency) RED, GREEN, YELLOW, VIOLET LED , 3D Vibration (Optional)	5 minutes	Daily	EMS Mode Apply skincare product to the skin and gently glide the device across the treatment area. Helps support facial contours and improve the appearance of skin firmness.

Intensity Selection Guide

- Almost no heat sensation ➔ **Normal**
- Warm sensation ➔ **Appropriate level**
- Feels hot ➔ **Lower the level**

* Results may vary depending on individual skin condition and usage environment.

Expected Experience Guide

- **After 1 use:** Skin may feel smoother and look more radiant
- **After 2–3 weeks:** Skin texture may appear more refined
- **After 4 weeks or more:** Skin may appear firmer and facial contours more defined

careCL Aesthetic Programs

careCL TECH FIT can be used in various skin care routines by combining the SHOT, GLOW, and EMS modes.

Mode
Automatic Time

SHOT Mode
approx.
7 minutes auto stop

GLOW Mode
approx.
5 minutes auto stop

EMS Mode
approx.
5 minutes auto stop

* Each mode automatically stops after the preset treatment time.

1. Skin Firmness Program

Cheeks, Nasolabial Area, Jawline

**FACE
SHOT**
7min.

Skin firmness care using the 5-second stamping technique

GLOW
5min.

Supports skincare absorption

EMS
5min.

Helps support skin firmness as a finishing step

2. V-Line Program

Jawline, Chin Area

**FACE
SHOT**
7min.

Focused care around the jawline

EMS
5min.

Skin firmness care around the jawline

3. Hydration Program

Full Face

GLOW
5min.

Supports skincare absorption

EMS
5min.

Supports skincare absorption

4. Fine Line Program

Eye Area, Lip Area

**EYE
SHOT**
7min.

Gentle stamping care around the eyes and lip area

GLOW
5min.

Supports absorption of eye cream or serum

* When using around the eyes, a lower intensity level is recommended.

5. Neck Program

Neck

**BODY
SHOT**
7min.

Focused care on the neck area

GLOW
5min.

Supports skincare absorption

EMS
5min.

helps support skin firmness

6. Hand Firmness Program

Back of Hands

GLOW
5min.

Supports absorption of hand cream

EMS
5min.

helps support skin firmness

7. Arm Firmness Program

Upper Arms

**BODY
SHOT**
7min.

Focused care on the arm area

EMS
5min.


helps support skin firmness

5 How to Use the Device

1) Preparation Before Use

- ① **Check Charging Status:** Make sure the battery is sufficiently charged.
(The device may stop automatically during use if the battery level is low.)

- ② **Cleansing:** Cleanse your face thoroughly with a cleanser and apply basic skincare products (such as toner and moisturizer).

③ **Skin Preparation**

- When using a mask pack: Take out the careCL Grid Mask Pack and attach it to the face.
Remove air bubbles from the center of the face toward the outer edges and smooth it without wrinkles.
Follow the large grid areas during use, and repeat on smaller grid areas for more detailed care.

- Apply careCL Collagen Gel evenly to the treatment area.
(Use a sufficient amount.)


TIP


EMS mode is recommended to be used with collagen gel. GLOW mode works well with water-based skincare products. Products with high oil content are not recommended.

2) Power & Basic Operation

- ① **Power On:** Press and hold the power / mode button for over 2 seconds.
- ② **Mode Selection:** Press the power / mode button briefly to cycle through
FACE SHOT ⌚ EYE SHOT ⌚ BODY SHOT ⌚ GLOW ⌚ EMS


- ③ **Intensity:** Press the level button to choose from LEVEL 1-5. (LEVEL 1 recommended for first-time users)
*** Start from LEVEL 3 when using with a mask pack**

- ④ **Vibration (Optional):** Press and hold the level button for 2 seconds to toggle vibration ON / OFF.

Output level /
Vibration control button

- ⑤ **Power Off:** Press and hold the power / mode button for over 2 seconds.

Power / Mode
button

Battery Charging


• Make sure the battery is fully charged before use

- ✓ If LOW BATTERY appears when powered ON, charge immediately.
- ✓ Use the dedicated included cable. Use only 5V 2A certified adapters.
- * Using a fast charger for phones may damage the product.

Charging in Progress Output level icons blink by step (1–5 levels)

Charging Complete Output level icon blinks at level 5, turns off after 3 minutes

* Device functions do not work during charging.

6 How to Use Each Mode


3) SHOT Mode General Instructions (FACE / EYE / BODY)

5-Second Stamping Shot

SHOT mode is designed with a "5-second stamping shot" technique that focuses energy on the skin by maintaining contact with one area, rather than rubbing across the skin.

- ① **Operation Method** : Stamping (press and hold). Maintain contact for about 5 seconds, then the device will automatically pause with a "ding" sound. Lift the device and move to the next area.
- ② **Duration**: About 7 minutes per area (Mode automatically stops after 7 minutes / repeat 2-3 times if needed).

A FACE SHOT (For Face Use)

Mode Selection FACE SHOT ☉ Start from LEVEL 3 when using with a mask pack.
Start from LEVEL 1 when using with collagen gel.

First Area Begin at the center of the cheek (below the cheekbone). Place the device head vertically on a grid section or a 2-3 cm area.

Stamping When you hear the "ding" sound, maintain gentle pressure for about 5 seconds. When the "ding" sound (pause) occurs, lift the device straight up.


Move Move to the next grid section and repeat, covering the entire face evenly.

Sequence Cheeks ☉ Nasolabial Area ☉ Jawline ☉ Between the Brows ☉ Forehead

Caution Do not use directly on the eyelids or mucous membranes. Avoid the front of the neck (thyroid area).

TIP Avoid repeating the same area more than 10 consecutive times. Perform daily care within 10-20 minutes (mode auto-stops after 7 minutes; repeat 2-3 times if needed).

B EYE SHOT (For Eye Area)

Mode / Level EYE SHOT, start from **LEVEL 1**

Area Use along the orbital bone around the eye. Do not approach within 2cm of the lash line.

Method Same as FACE mode – stamp for 5 seconds, then move to the next area.

Time Alternate between both eyes within 7 minutes.

C BODY SHOT (For Body Use)


Mode / Level BODY SHOT, start from **LEVEL 1~3**(Adjust depending on area thickness/sensitivity).

Method Neck sides, forearms, backs of hands, abdomen, thighs, etc.

Time Use the same **5-second stamping** technique as in FACE mode to systematically care for all areas.

4) GLOW Mode (Glow & Texture Care)

- ❖ Purpose: Improves glow, radiance, and skin texture.

Application Apply a thin layer of water-based serum/gel (Do not use oil-based products).


Mode / Level Select GLOW, start from **LEVEL 1**.

Method Gently glide over entire face in the order: Forehead ☉ Cheeks ☉ Jawline.

Time 5–10 minutes

5) EMS Mode (Muscle Stimulation)

- ❖ Purpose: Activate muscles and improve firmness using microcurrents.

Apply Gel Apply a generous amount of collagen gel to the treatment area.

Mode / Level Select EMS mode, start from **LEVEL 1**

Contact When device contacts the skin, it will activate with a “beep” (no 5-second pause like SHOT mode).

Method Slow glide: Move slowly along the direction of muscle fibers (e.g., jawline ☉ ear).

Time Use for 10 minutes in total

Caution Do not use on those with heart conditions, metal plates, or implants. Avoid front neck, wounded, or inflamed areas. If tingling is uncomfortable, lower the level or stop use.

6) Finishing & Aftercare

Power Off Press and hold the Power/Mode button for about 2 seconds to turn the device OFF. After use, make sure the device is completely powered off.

Device Care After use, gently wipe off any remaining gel or essence on the device head with a soft dry cloth. Contact with water may cause product damage. For long-term storage, turn off the device completely before storing. Handle the device carefully to avoid impact.

Skincare After using SHOT or EMS mode, you can use GLOW mode to help gently absorb any remaining gel or essence on the skin. Avoid strong facial massage, exfoliation, or stimulation immediately after use. If your skin feels dry, apply sufficient moisture using a toner or cream. When going outdoors, it is recommended to apply sunscreen.

Product Specifications

Product	careCL TECH FIT BEAUTY DEVICE
Included	1 Device, 1 Cradle, 1 USB Cable
Material	PC, ABS, Stainless Steel
Size	Approx. 41(W) x 40(D) x 178(H) mm
Weight	Approx. 180g
Power	DC 3.7V 3300mAh Li-ion Rechargeable Battery Usage Time: 150 minutes, Charging Time: 180 minutes
Patents & Certifications	Technology patent, Design patent    certified
Country of Manufacture	South Korea
Manufacturer	Lamu Inc.

7 Storage Method & Precautions

How to Store

- After each use, wipe off any product residue with a soft cloth.
- Avoid dropping the device.
- Store away from direct sunlight.
- Avoid storing in places with chemicals, gas, or hazardous substances.
- Do not store in places below 10°C or above 50°C as this may cause malfunction.
- Keep away from corrosive gases or humidity.

* The manufacturer is not responsible for damage or issues caused by user negligence or misuse.

Precautions

⊗ Major Warnings (Do not use if you are)

- A user with a pacemaker, defibrillator, or implanted electronic device
- A child / pregnant woman
- Someone with open wounds or hypertrophic scars
- Someone with severe acne or extremely sensitive skin
- Wearing any implanted medical device
- Having systemic or local skin conditions that affect wound healing

✓ Cautions During Use

- Read the user guide carefully before use.
- Adjust intensity based on skin sensitivity.
- Always apply sufficient careCL Collagen Gel or use the mask pack properly attached to the skin to ensure full contact.
(Device will not operate unless in proper contact with skin.)
- Do not use over areas where filler or autologous fat was recently injected (wait a set period before use).
- Do not use directly on eyeballs or upper eyelids as it may cause serious damage including blindness.
- Read the prohibited usage items carefully. If any apply, do not use.

✓ Cautions After Use

- Redness or warmth may appear on the treated facial area after use.
(This is a normal reaction and usually disappears within a few hours.)
- Avoid saunas after use.
- Do not irritate the skin post-treatment.
- Even if the skin itches, avoid scratching or touching the treated area.

⊗ Forbidden Areas for careCL TECH FIT

- Do not use on eyes or lips directly. Avoid upper eyelids, Adam's apple, etc.
- Always apply sufficient careCL Collagen Gel or attach the mask pack firmly for full contact with the skin.

Frequently Asked Questions (FAQ)

Q Why should the device be used for 5 seconds on one area?

careCL TECH FIT uses a stamping method, where the device is held on one spot for a short period of time. Home-use devices are designed within safety standards. Maintaining contact on single area allows the device to deliver energy in a stable and consistent manner. Through internal testing, 5 seconds was determined to be an appropriate duration that balances effective energy delivery with user comfort. Rather than quickly moving across the skin, briefly focusing on one area at a time is the key concept of the TECH FIT usage method.

Q What is the role of the Grid Mask Pack?

The Grid Mask Pack serves as a visual guide that helps users easily identify the recommended areas for device use. When used on its own, it also functions as a hydrating and conditioning mask pack designed to support overall skin care.

Key ingredients include	• Peptide complex • Multi-hyaluronic acid complex • Galactomyces ferment filtrate These ingredients help support skin hydration and overall skin condition.
-------------------------	--

Q Are there any additional consumable costs?

All modes of careCL TECH FIT can be used without additional consumables or replacement cartridges. The device is designed so that all modes can be used repeatedly without additional running costs.

Q Why does the device not operate when it is not touching the skin?

The device is equipped with a skin contact sensor. Energy output is activated only when proper skin contact is detected, which helps improve safety during use.

Q What is the difference between SHOT / GLOW / EMS modes?	
SHOT Mode (Stamping Energy Technology)	• FACE – Focused care for facial areas (RED LED) • EYE – Focused care for the eye area (GREEN LED) • BODY – Care for neck and body areas (YELLOW LED) Different settings and LED indicators are applied depending on the treatment area.
GLOW Mode	GLOW mode uses EP (electroporation) / galvanic technology. When used after applying skincare products, it can help support the absorption of cosmetic products and improve skin condition.
EMS Mode	EMS mode uses microcurrent stimulation and can help support muscle stimulation and skin care in areas such as the jawline and facial contours.

Q Which level should I start with when using the device for the first time?

For sensitive skin, it is recommended to start with LEVEL 1. Once your skin becomes accustomed to the device, you can gradually increase the level as needed.

Q Is it normal to feel warmth while using the device?

When the careCL **TECH FIT** device operates, you may feel a gentle warming sensation on the skin. This is a normal characteristic of the device.
The level of warmth felt may vary depending on individual skin condition and sensitivity.

Q How should I store the device after use?

If cosmetic products remain on the device, wipe the surface with a soft cloth or tissue.
After cleaning, store the device in a dry place.

Q Can I use the device while charging?

For safety reasons, the device does not operate while charging.
Please use the device after charging is complete.

Q Is it normal to feel a tingling sensation in EMS mode?

Depending on the selected level and skin condition, a mild tingling sensation may be felt during EMS mode.
If the sensation feels uncomfortable, reduce the level or stop using the device.

Q Can the voice alert or sound function be turned off?

The voice alert and sound notification are basic features provided by the device specifications and cannot be turned off.

Q Can the device be used during pregnancy?

Use during pregnancy or breastfeeding is not recommended.
Due to hormonal changes, temporary skin sensitivity or skin concerns may occur.

Q What is the battery operating time?

Battery usage time may vary depending on the selected mode and level.
When the battery is low, the LED indicator will blink red five times. Please recharge the device.

Q How do I turn the vibration function on or off?

Press and hold the level button for approximately 2 seconds to turn the vibration function ON or OFF.

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